



Ready-to-use 3D Human Mini Liver Spheroids (SP3D-HMLS)

Catalog #SP3D-5108

Product Description

The liver is a vital metabolic organ and its functional unit is termed hepatic lobule, which is comprised of hexagonally-arranged hepatocytes innervated by the hepatic sinusoids. The liver has a multi-dimensional role, including protein synthesis, nutrient storage, bile production and bilirubin metabolism. Another key hepatic function is the clearance of first pass doses of orally-absorbed substances, such as drugs and xenobiotics. Drug-induced hepatotoxicity remains a major concern for pharmaceutical industries. Presently, 2D culture of primary hepatocytes is the gold standard for hepatotoxicity testing. This approach, however, is limited by factors such as poor proliferative capacity of primary hepatocytes and their tendency to rapidly dedifferentiate in 2D culture [1]. Furthermore, the cellular diversity and liver architecture are not taken into consideration. The crosstalk between hepatocytes and the non-parenchymal cells in the liver, such as sinusoidal endothelial cells and hepatic stellate cells, significantly influences hepatic function and the accuracy of predictive hepatotoxicity testing.

To provide an *in vitro* liver model that better reflects hepatic function physiologically, ScienCell has developed ready-to-use human mini liver spheroids comprised of human hepatocytes, hepatic stellate cells, and sinusoidal endothelial cells. The co-culture of these main hepatic cell types allows the functional maintenance of cells over time due to a more native, three-dimensional (3D) architecture. These spheroids are ready for experiments within 24-48 hours of thawing, and are an excellent *in vitro* model for studying liver function and for the examination of hepatotoxicity.

Kit Components (Included)

3D Cell Culture Components				
Cat #	# of vials	Product Name	Quantity	Storage
SP-5108	1	Human Mini Liver Spheroids (SP-HMLS)	4 × 10 ³ spheroids	Liquid nitrogen
3D-5201	1	3D-Liver Spheroid Medium (3D-LSpM)	200 mL	2-8 °C
3D-5252	1	3D-Liver Spheroid Supplement (3D-LSpS)	4 mL	-20 °C
0583	1	Penicillin/Streptomycin Solution (P/S)	2 mL	-20 °C
0343 (or) 0353 (or) 0383	1	Ultra-Low Binding Culture Plates (24-, 48-, or 96- well plate)	1 plate	RT

Quality Control

SP3D-HMLS are tested for the formation of functional and uniform 3D human mini liver spheroids according to the included protocol. All components are negative for bacterial and fungal contamination.

Rev. 0

Product Use

SP3D-HMLS are for research use only. It is not approved for human or animal use, or for application in clinical or *in vitro* diagnostic procedures.

Shipping

SP-5108, 3D-5252, 0583 are shipped on dry ice. 3D-5201 and (0343 or 0353 or 0383) are shipped at room temperature.

References

[1] Elaut G, Henkens T, Papeleu P, Snykers S, Vinken M, Vanhaecke T, Rogiers V. (2006) "Molecular mechanisms underlying the dedifferentiation process of isolated hepatocytes and their cultures." *Curr Drug Metab.* 7: 629-660.

Procedure:

Step I: Preparing the complete 3D culture medium

1. Thaw 3D-liver spheroid supplement (3D-LSpS; Cat. #3D-5252), and penicillin/streptomycin solution (P/S solution; Cat. #0583) at 37°C. Mix 3D-LSpS and P/S solution into the 3D-liver spheroid medium (3D-LSpM; Cat. #3D-5201) by gently swirling the medium bottle around.
 - a. 3D-LSpM is **viscous** and optimized for homogenous spheroid formation.
 - b. Warm the complete 3D-LSpM to **room temperature** before use.
 - c. When stored in the dark at 4°C, the complete medium is stable for one month.

Step II: Thawing and maintaining the ready-to-use 3D spheroids

2. One frozen vial contains $\geq 4 \times 10^3$ spheroids, which is sufficient for plating into half of a multi-well plate (e.g. 24-, 48-, and 96-well ultra-low binding culture plate).
3. Place the frozen vial in a 37°C water bath. Hold and rotate the vial gently until the contents completely thaw. Promptly remove the vial from the water bath, wipe it down with 70% ethanol, and transfer it to the sterile field.
4. Carefully remove the cap without touching the interior threads. Gently pipette the spheroid suspension up and down **two times** to disperse potential spheroid aggregates.
5. Gently transfer the spheroid suspension into a fresh 50 mL conical tube.
6. Add 12 mL of 3D culture medium to the above 50 mL conical tube.
7. Resuspend spheroids in 3D culture medium by gently pipetting up and down for ~ 5 times using a serological pipette.

Note: 3D culture medium has a high viscosity; thus, pipetting slowly is important to avoid bubble formation.

