



**Small Airway Epithelial Cell Medium-basal-phenol red free
(SAEpiCM-b-prf)
Catalog #3231-b-prf**

Product Description

Small Airway Epithelial Cell Medium-basal-phenol red free (SAEpiCM-b-prf), when used with Small Airway Epithelial Cell Growth Supplement (SAEpiCSG, Cat #3272), is a complete medium designed for optimal growth of normal small airway epithelial cells *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds and trace minerals. The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO₂/95% air.

Components

SAEpiCM-b-prf consists of 500 ml of basal medium.

Product Use

SAEpiCM-b-prf is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

Storage

Store at 4°C. Protect from light.

Shipping

Room temperature.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.