



**Human Esophageal Smooth Muscle Cell Pellet  
(HEsSMC Cell Pellet)**  
Catalog #CP2710

Quantity: 5 million cells

**Product Description**

Human Esophageal Smooth Muscle Cell Pellet (HEsSMC Cell Pellet) is prepared from early passage Human Esophageal Smooth Muscle Cells. Each pellet contains 5 million cells and can be used for a variety of applications including PCR, western blotting, genomic DNA library construction, and gene expression profiling. Cell pellets from ScienCell Research Laboratories are convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

**Applications**

- Genomic DNA
- cDNA
- Total RNA
- microRNA
- Protein Lysates

**Product Use**

HEsSMC Cell Pellet is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

**Storage**

Store cell pellet immediately at -80°C upon arrival. Avoid freeze/thaw cycles.

**Shipping**

Dry ice.