



Buffer VL
(BuffVL)
Catalog #MB891a
35 mL

Description

ScienCell's Buffer VL is a component of Viral RNA Extraction Kit (MB891), which provides a fast and reliable way to purify viral RNA from a variety of samples such as cell culture media, plasma, serum, saliva, nasal/nasopharyngeal/oropharyngeal swab specimens, bronchoalveolar lavage, and other body fluids. The Buffer VL is optimized for the Viral RNA Extraction Kit.

Quality control

BuffVL is tested as a component of Viral RNA Extraction Kit (MB891).

Product Use

BuffVL is for research use only. It is not approved for human or animal use, or for application in clinical or *in vitro* diagnostic procedures. BuffVL is only optimized for Viral RNA Extraction Kit (MB891). Please refer the product sheet of Viral RNA Extraction Kit (MB891) for detailed procedures.

Shipping and Storage

Ambient temperature.